## gmm ver2

June 7, 2025

```
[73]: import numpy as np
      import torch
      import matplotlib.pyplot as plt
      import scipy.signal as signal
      import torch.distributions as dist
      from jaxtyping import Float, Int
      from torch import Tensor
      import spikeinterface.extractors as se
      import spikeinterface.comparison as sc
      from spikeinterface.preprocessing import bandpass_filter
      from sklearn.metrics import confusion_matrix, adjusted_rand_score
      from scipy.signal import find_peaks
      import h5py
[17]: def load_mearec_data(recording_file="recordings.h5"):
          Load MEArec data and convert to PyTorch tensors compatible with the \Box
       ⇔notebook format
          11 11 11
          # Load MEArec data
          recording, sorting_true = se.read_mearec(recording_file)
          # Apply preprocessing (bandpass filter)
          recording_filt = bandpass_filter(recording, freq_min=300, freq_max=6000)
          # Get traces as numpy array and convert to torch tensor
          traces = torch.tensor(recording_filt.get_traces(), dtype=torch.float32)
          # Get ground truth spike times and IDs
          true spike times = []
          true_spike_ids = []
          for unit_id in sorting_true.unit_ids:
              unit_spikes = sorting_true.get_unit_spike_train(unit_id)
              true_spike_times.extend(unit_spikes)
              # Convert unit id to integer (remove '#' if present)
```

```
unit_idx = int(unit_id.lstrip('#')) if isinstance(unit_id, str) else__

    int(unit_id)

             true_spike_ids.extend([unit_idx] * len(unit_spikes))
         # Sort by spike times
         sorted indices = np.argsort(true spike times)
         true_spike_times = torch.tensor([true_spike_times[i] for i in_
       ⇔sorted_indices], dtype=torch.long)
         true_spike_ids = torch.tensor([true_spike_ids[i] for i in sorted_indices],_
       →dtype=torch.long)
         # Extract true templates from MEArec data
         with h5py.File(recording_file, 'r') as f:
             templates = f['templates'][:] # Shape: (n_units, n_jitters, __
       \rightarrow n channels, n samples)
         # Convert templates to the expected format: (n_units, n_samples, n_channels)
         # Take median across jitters and transpose
         →1), dtype=torch.float32)
         # Constants
         sampling_frequency = recording.sampling_frequency
         num_samples, num_channels = traces.shape
         true_num_spikes = len(true_spike_times)
         true_num_neurons, window_size, _ = true_templates.shape
         return {
             'traces': traces,
             'true_spike_times': true_spike_times,
             'true_spike_ids': true_spike_ids,
             'true_templates': true_templates,
             'sampling_frequency': sampling_frequency,
             'num_samples': num_samples,
             'num_channels': num_channels,
             'true_num_spikes': true_num_spikes,
             'true_num_neurons': true_num_neurons,
             'window_size': window_size
         }
[98]: data = load_mearec_data("data/sim_data/recordings5.h5")
      # Extract variables in the same format as the original notebook
     traces = data['traces']
     true_spike_times = data['true_spike_times']
     true_spike_ids = data['true_spike_ids']
```

true\_templates = data['true\_templates']

```
SAMPLING_FREQUENCY = data['sampling_frequency']
      NUM_SAMPLES, NUM_CHANNELS = traces.shape
      TRUE_NUM_SPIKES = data['true_num_spikes']
      TRUE_NUM_NEURONS, WINDOW_SIZE, _ = true_templates.shape
      WINDOW_SIZE = 96
      print(f" Loaded MEArec data:")
      print(f" Traces shape: {traces.shape}")
      print(f" True spikes: {TRUE_NUM_SPIKES}")
      print(f" True neurons: {TRUE NUM NEURONS}")
      print(f" Sampling frequency: {SAMPLING_FREQUENCY} Hz")
      print(f" Window size: {WINDOW SIZE}")
      Loaded MEArec data:
        Traces shape: torch.Size([640000, 4])
        True spikes: 669
        True neurons: 4
        Sampling frequency: 32000.0 Hz
        Window size: 96
[78]: true_templates.shape
[78]: torch.Size([8, 416, 32])
[79]: def plot_slice(traces: Float[Tensor, "num_samples num_channels"],
                     spike_times: Int[Tensor, "num_spikes"],
                     spike_ids: Int[Tensor, "num_spikes"],
                     start_time: float=0,
                     end time: float=1.0,
                     ylim: float=None):
          11 11 11
          Plot a slice of the recording and the corresponding spike times.
          # Extract constants
          num_samples, num_channels = traces.shape
          num_neurons = spike_ids.max() + 1
          # Extract the slice of the recording
          slc = slice(int(start_time * SAMPLING_FREQUENCY), int(end_time *_
       →SAMPLING_FREQUENCY))
          # Get the y-limit for traces if not specified
          if ylim is None:
              ylim = 1.05 * traces[slc].abs().max()
          # Make an array of times in seconds
          times_sec = torch.arange(0, num_samples) / SAMPLING_FREQUENCY
```

```
spike_in_slice = (spike_times_sec >= start_time) & (spike_times_sec <=_
       ⇔end_time)
         # Plot the spike raster
         fig, axs = plt.subplots(num channels + 1, 1, sharex=True, figsize=(6, 1.5 *
       →(num channels + 1)))
         axs[0].plot(spike_times_sec[spike_in_slice], spike_ids[spike_in_slice],_u
       axs[0].set_ylabel("neuron")
         axs[0].set_yticks(torch.arange(0, num_neurons))
         axs[0].set_title("spike train")
         # Plot the voltage traces
         for i in range(num_channels):
             axs[i+1].plot(times_sec[slc], traces[slc, i], color='k', lw=1, alpha=0.
       ⇒5)
             axs[i+1].set_title(f"channel {i+1}", fontsize=10)
             axs[i+1].set ylim(-ylim, ylim)
             axs[i+1].set_ylabel("trace\n[z-score]")
             axs[i+1].set xlim(start time, end time)
             if i == num channels - 1:
                 axs[i+1].set_xlabel(f"time [sec]")
         plt.tight_layout()
         return fig, axs
[80]: def plot_waveforms(waveforms: Float[Tensor, "num_waveforms windows_size_
       name: str="waveform"):
         Plot a grid of spike waveforms or templates.
         num_waveforms, window_size, num_channels = waveforms.shape
         times = torch.arange(-window_size//2, window_size//2) / SAMPLING_FREQUENCY_
       →* 1000
         fig, axs = plt.subplots(num_channels, num_waveforms,
                                figsize=(1.25 * num_waveforms, 1.25 * num_channels),
                                sharex=True, sharey=True)
          # Handle case where there's only one channel or one waveform
         if num_channels == 1:
             axs = axs.reshape(1, -1)
         if num_waveforms == 1:
             axs = axs.reshape(-1, 1)
```

spike\_times\_sec = times\_sec[spike\_times]

```
lim = 1.05 * waveforms.abs().max()
          for j in range(num_waveforms):
              for i in range(num_channels):
                  axs[i,j].plot(times, waveforms[j,:,i], color='k', lw=1, alpha=0.5)
                  axs[i,j].set_ylim(-lim, lim)
                  if j == 0:
                      axs[i,j].set_ylabel(f"channel {i}\n [z-score]")
                  if i == 0:
                      axs[i,j].set_title(f"{name} {j}")
                  if i == num channels - 1:
                      axs[i,j].set_xlabel(f"time [ms]")
          plt.tight_layout()
          return fig, axs
[81]: def evaluate_detected_spikes(
          detected_spike_times: Int[Tensor, "num_detected_spikes"],
          true_spike_times: Int[Tensor, "num_true_spikes"]) -> dict:
          Evaluate detected spikes against ground truth with boundary handling
          num_windows = NUM_SAMPLES // WINDOW_SIZE
          # Filter out spikes that would cause out-of-bounds errors
          # Keep only spikes that fall within valid window boundaries
          valid_true_spikes = true_spike_times[true_spike_times < num_windows *_u
       →WINDOW_SIZE]
          valid_detected_spikes = detected_spike_times[detected_spike_times <_
       →num_windows * WINDOW_SIZE]
          # Create masks for windows containing spikes
          true_spike_mask = torch.zeros(num_windows, dtype=torch.bool)
          if len(valid_true_spikes) > 0:
              true_spike_mask[valid_true_spikes // WINDOW_SIZE] = True
          detected_spike mask = torch.zeros(num_windows, dtype=torch.bool)
          if len(valid_detected_spikes) > 0:
              detected_spike_mask[valid_detected_spikes // WINDOW_SIZE] = True
          # Calculate confusion matrix elements
          tp = (true_spike_mask & detected_spike_mask).sum().item()
          fp = (detected_spike_mask & ~true_spike_mask).sum().item()
          fn = (~detected_spike_mask & true_spike_mask).sum().item()
```

tn = (~detected\_spike\_mask & ~true\_spike\_mask).sum().item()


```
return dict(tp=tp, fp=fp, fn=fn, tn=tn)
```

```
[82]: class SimpleSpikeSorter:
          A simple spike sorter that uses a Gaussian mixture model (GMM) to fit the \Box
       ⇔templates.
          HHHH
          def __init__(self, initial_templates: Float[Tensor, "num_neurons⊔
       ⇔window_size num_channels"],
                       noise scale: float = 1.0) -> None:
              self.templates = initial_templates
              self.noise_scale = noise_scale
          @property
          def num_neurons(self) -> int:
              return self.templates.shape[0]
          def log probability(self, spike_waveforms: Float[Tensor, "num_spikes_
       ⇔window_size num_channels"],
                             spike_ids: Int[Tensor, "num_spikes"]) -> float:
              """Compute the log probability of the spike waveforms given the

→templates"""
              assigned_templates = self.templates[spike_ids]
              11 = torch.distributions.Normal(assigned_templates, self.noise_scale).
       →log_prob(spike_waveforms)
              11 = 11.sum().item()
              return 11
          def update_templates(self, spike_waveforms: Float[Tensor, "num_spikes_
       ⇔window size num channels"],
                              spike_ids: Int[Tensor, "num_spikes"]):
              """Update the templates using the spike waveforms and spike IDs"""
              for k in range(self.num_neurons):
                  template_spikes = spike_waveforms[spike_ids == k]
                  if len(template_spikes) > 0:
                      self.templates[k] = template_spikes.mean(dim=0)
          def update_spike_ids(self, spike_waveforms: Float[Tensor, "num_spikes_u
       ⇔window_size num_channels"]
                              ) -> Int[Tensor, "num_spikes"]:
              """Update the spike IDs using the templates and spike waveforms"""
              template_reshaped = self.templates.unsqueeze(1)
              waveforms_reshaped = spike_waveforms.unsqueeze(0)
              11 = torch.distributions.Normal(template reshaped, self.noise scale).
       →log_prob(waveforms_reshaped)
```


```
spike_ids = 11.sum(dim=(2,3)).argmax(dim=0)
      return spike_ids
  def fit(self, spike_waveforms: Float[Tensor, "num_spikes_window_size_
num_iterations: int=10) -> tuple[Float[Tensor, "num_iterations"],__
→Int[Tensor, "num_spikes"]]:
      """Fit the GMM to the spike waveforms using coordinate ascent"""
      # Initialize the spike IDs
      spike_ids = self.update_spike_ids(spike_waveforms)
      # Iterate to update the templates and spike IDs
      lps = []
      for itr in range(num_iterations):
          lps.append(self.log_probability(spike_waveforms, spike_ids))
          print(f"iteration {itr}: log probability = {lps[-1]:.2f}")
          # Update the templates and spike IDs
          self.update_templates(spike_waveforms, spike_ids)
          spike_ids = self.update_spike_ids(spike_waveforms)
      return torch.tensor(lps), spike_ids
```

```
[99]: # Plot the first second of MEArec data
_ = plot_slice(traces, true_spike_times, true_spike_ids, start_time=0.0, u
→end_time=1.0)
```



```
[100]: # Plot true templates
fig, axs = plot_waveforms(true_templates, name="template")
fig.suptitle("MEArec True Templates", y=1.05)
plt.show()
```

## **MEArec True Templates**



```
# Convert back to torch tensors
si_spike_times = torch.tensor(peaks['sample_index'])
print(f" SpikeInterface detected {len(si_spike_times)} spikes")

# Extract waveforms using SpikeInterface detection
#spike_waveforms_si = torch.zeros((len(si_spike_times), WINDOW_SIZE, USINDOW_CHANNELS))
traces_np = recording_filt.get_traces()

# for i, t in enumerate(si_spike_times):
# start_idx = t - WINDOW_SIZE // 2
# end_idx = t + WINDOW_SIZE // 2

# if start_idx >= 0 and end_idx < traces_np.shape[0]:
# spike_waveforms_si[i] = torch.tensor(traces_np[start_idx:end_idx])</pre>
```

Using SpikeInterface peak detection...

```
noise_level (no parallelization): 0%| | 0/20 [00:00<?, ?it/s] detect peaks using locally_exclusive (no parallelization): 0%| | 0/20_ _{\hookrightarrow} [00:00<?, ?it/s]
```


SpikeInterface detected 653 spikes

```
[103]: def extract_waveforms_robust(traces_np, spike_times, window_size):
           Robust waveform extraction with proper boundary handling
           num_spikes = len(spike_times)
           num_samples, num_channels = traces_np.shape
           half_window = window_size // 2
           spike_waveforms = []
           valid_spike_times = []
           print(f"Extracting waveforms with window size {window size}...")
           for i, t in enumerate(spike_times):
               start_idx = int(t - half_window)
               end_idx = int(t + half_window)
               # Check boundaries
               if start_idx >= 0 and end_idx < num_samples:</pre>
                   waveform = traces_np[start_idx:end_idx, :]
                   # Verify we got the right shape
                   if waveform.shape[0] == window_size:
```

```
spike_waveforms.append(torch.tensor(waveform))
                valid_spike_times.append(t)
            else:
                print(f"Skipping spike {i}: wrong shape {waveform.shape}")
        else:
            print(f"Skipping spike {i}: out of bounds (start={start_idx},__
 →end={end_idx})")
   if len(spike_waveforms) > 0:
        spike_waveforms = torch.stack(spike_waveforms)
       print(f" Extracted {len(spike_waveforms)} valid waveforms")
   else:
       print(" No valid waveforms extracted!")
        spike_waveforms = torch.zeros((0, window_size, num_channels))
   return spike_waveforms, torch.tensor(valid_spike_times)
# Apply robust extraction
spike_waveforms_fixed, valid_spike_times = extract_waveforms_robust(
   traces_np, si_spike_times, WINDOW_SIZE
)
# Plot the results
if len(spike_waveforms_fixed) > 0:
   fig, axs = plot_waveforms(spike_waveforms_fixed[:5], name="waveform")
   fig.suptitle("Robust Extraction - First 5 Waveforms", y=1.05)
   plt.show()
else:
   print("No waveforms to plot!")
```

Extracting waveforms with window size 96... Extracted 653 valid waveforms


## Robust Extraction - First 5 Waveforms




```
plt.figure(figsize=(8, 4))
plt.plot(lps, '-o')
plt.xlabel("iteration")
plt.ylabel("log probability")
plt.title("GMM Convergence")
plt.grid(True, alpha=0.3)
plt.show()

# Plot inferred templates
fig, axs = plot_waveforms(spike_sorter.templates, name="template")
fig.suptitle("GMM Inferred Templates", y=1.05)
plt.show()
```

Running GMM spike sorting with 4 neurons...
iteration 0: log probability = -3749877.25
iteration 1: log probability = -2579594.50
iteration 2: log probability = -2501290.75
iteration 3: log probability = -2420767.00
iteration 4: log probability = -2415596.75
iteration 5: log probability = -2415596.75
iteration 6: log probability = -2415596.75
iteration 7: log probability = -2415596.75
iteration 8: log probability = -2415596.75
iteration 9: log probability = -2415596.75



## **GMM Inferred Templates**




```
[105]: # Alternative: Use ground truth spike times for cleaner evaluation
def evaluate_gmm_with_ground_truth(true_spike_times, true_spike_ids, traces,
gmm_sorter):
    """
        Evaluate GMM using ground truth spike times
    """
        print(" Evaluating GMM using ground truth spike times...")

# Extract waveforms at ground truth locations
gt_waveforms = []
valid_gt_times = []
valid_gt_times = []
for i, t in enumerate(true_spike_times):
        start_idx = t - WINDOW_SIZE // 2
```

```
end_idx = t + WINDOW_SIZE // 2
        if start_idx >= 0 and end_idx < traces.shape[0]:</pre>
            waveform = traces[start_idx:end_idx]
            gt_waveforms.append(waveform)
            valid_gt_times.append(t)
            valid_gt_ids.append(true_spike_ids[i])
    gt waveforms = torch.stack(gt waveforms)
    valid_gt_ids = torch.tensor(valid_gt_ids)
    print(f" Extracted {len(gt_waveforms)} ground truth waveforms")
    # Apply GMM to ground truth waveforms
    gmm_predictions = gmm_sorter.update_spike_ids(gt_waveforms)
    # Create confusion matrix
    cm_gt = confusion_matrix(valid_gt_ids.numpy(), gmm_predictions.numpy())
    return cm_gt, valid_gt_ids, gmm_predictions
# Apply ground truth evaluation
cm_gt, gt_ids, gmm_pred = evaluate_gmm_with_ground_truth(
    true_spike_times, true_spike_ids, traces, spike_sorter
)
# Plot ground truth-based confusion matrix
plt.figure(figsize=(10, 8))
plt.imshow(cm_gt, interpolation="none", cmap=plt.cm.Blues)
plt.xticks(range(cm_gt.shape[1]))
plt.yticks(range(cm_gt.shape[0]))
plt.colorbar()
plt.xlabel('GMM Predicted Clusters')
plt.ylabel('True Neuron IDs')
plt.title("GMM Performance on Ground Truth Spike Times")
for i in range(cm_gt.shape[0]):
    for j in range(cm_gt.shape[1]):
        plt.text(j, i, str(cm_gt[i, j]), ha='center', va='center',
                color='white' if cm_gt[i, j] > cm_gt.max() / 2 else 'black')
plt.tight_layout()
plt.show()
# Calculate ARI for ground truth evaluation
ari_gt = adjusted_rand_score(gt_ids.numpy(), gmm_pred.numpy())
print(f" GMM Performance on Ground Truth Spikes:")
```

```
print(f" Adjusted Rand Index: {ari_gt:.3f}")
```

Evaluating GMM using ground truth spike times... Extracted 669 ground truth waveforms



GMM Performance on Ground Truth Spikes: Adjusted Rand Index: 0.973

```
[38]: import numpy as np
import torch
from fastdtw import fastdtw
from sklearn.decomposition import PCA
from sklearn.preprocessing import StandardScaler
from sklearn.mixture import GaussianMixture
from sklearn.metrics import silhouette_score
from joblib import Parallel, delayed
import time
```

```
from tqdm import tqdm
class DTWGMMSpikeSorter:
    HHHH
    DTW-based feature extraction with GMM clustering for spike sorting
    Based on:
    - PMC4749467: Dynamic Time Warping for spike classification
    - Nature Scientific Reports 2019: GMM overclustering approach
    - Hybrid methodology: DTW features + GMM clustering
    n n n
    def __init__(self, window_size: int, sampling_frequency: float,
                 n_reference_spikes: int = 50, dtw_radius: int = 3,
                 n_{jobs}: int = 4):
        self.window_size = window_size
        self.sampling_frequency = sampling_frequency
        self.n_reference_spikes = n_reference_spikes
        self.dtw_radius = dtw_radius
        self.n_jobs = n_jobs
        self.reference_spikes = None
        self.scaler = StandardScaler()
        self.pca = PCA(n_components=0.95) # Retain 95% variance
        self.gmm = None
    def _select_reference_spikes(self, spike_waveforms):
        Select diverse reference spikes using k-means++ initialization strategy
        From PMC4749467: "template estimator first calculates the averaged_{\sqcup}
 ⇒spike waveform"
        11 11 11
        n_spikes = len(spike_waveforms)
        if n_spikes <= self.n_reference_spikes:</pre>
            return spike waveforms
        # Flatten waveforms for distance computation
        flat_waveforms = spike_waveforms.reshape(n_spikes, -1)
        # K-means++ style selection for diverse references
        reference_indices = []
        # First reference: random selection
        reference_indices.append(np.random.randint(0, n_spikes))
        # Subsequent references: maximize minimum distance
        for _ in range(self.n_reference_spikes - 1):
            distances = []
```

```
for i in range(n_spikes):
              if i in reference_indices:
                   distances.append(0)
                   continue
              # Compute minimum distance to existing references
              min dist = float('inf')
              for ref_idx in reference_indices:
                   dist = np.linalg.norm(flat waveforms[i] -___
⇔flat_waveforms[ref_idx])
                   min_dist = min(min_dist, dist)
              distances.append(min_dist)
          # Select spike with maximum minimum distance
          next_ref = np.argmax(distances)
          reference_indices.append(next_ref)
      return spike_waveforms[reference_indices]
  def _compute_dtw_features_parallel(self, spike_waveforms):
      Compute DTW-based feature vectors using parallel processing
      From PMC4749467: "fastDTW method with automatic thresholding"
      n_spikes = len(spike_waveforms)
      n_refs = len(self.reference_spikes)
      print(f" Computing DTW features: {n_spikes} spikes × {n_refs}_{L}

→references...")
      def compute_dtw_row(spike_idx):
           """Compute DTW distances for one spike against all references"""
          spike = spike_waveforms[spike_idx]
          dtw row = np.zeros(n refs)
          for ref_idx, ref_spike in enumerate(self.reference_spikes):
               # Normalize spikes for DTW (from PMC4749467)
              spike_norm = spike / (np.linalg.norm(spike) + 1e-8)
              ref_norm = ref_spike / (np.linalg.norm(ref_spike) + 1e-8)
              # Compute fastDTW with radius constraint
              distance, _ = fastdtw(
                   spike_norm.flatten(),
                   ref_norm.flatten(),
                  radius=self.dtw_radius,
                  dist=lambda x, y: (x - y) ** 2
```

```
dtw_row[ref_idx] = distance
          return dtw_row
      # Parallel DTW computation
      dtw_features = Parallel(n_jobs=self.n_jobs, backend='threading')(
          delayed(compute_dtw_row)(i) for i in tqdm(range(n_spikes),__

desc="DTW features")

      return np.array(dtw_features)
  def _estimate_clusters_gmm_overclustering(self, features):
      Estimate optimal number of clusters using GMM overclustering approach
      From Nature Scientific Reports 2019: "initial overclustering of the \Box
⇔data"
      n_samples = len(features)
      max_clusters = min(20, n_samples // 10) # Reasonable upper bound
      print(f" Estimating clusters using GMM overclustering (max:

√{max_clusters})...")
      # Test different numbers of components
      bic_scores = []
      aic_scores = []
      silhouette_scores = []
      cluster_range = range(2, max_clusters + 1)
      for n_components in tqdm(cluster_range, desc="Testing cluster counts"):
           # Fit GMM with current number of components
          gmm = GaussianMixture(
              n_components=n_components,
              covariance_type='full',
              random_state=42,
              max iter=100
          )
          try:
              gmm.fit(features)
              labels = gmm.predict(features)
              bic_scores.append(gmm.bic(features))
               aic_scores.append(gmm.aic(features))
```

```
# Silhouette score for cluster quality
            if len(np.unique(labels)) > 1:
                sil_score = silhouette_score(features, labels)
                silhouette_scores.append(sil_score)
            else:
                silhouette_scores.append(-1)
        except Exception as e:
            print(f"Failed for {n_components} components: {e}")
            bic_scores.append(float('inf'))
            aic scores.append(float('inf'))
            silhouette_scores.append(-1)
    # Find optimal number using multiple criteria
   bic_optimal = cluster_range[np.argmin(bic_scores)]
    aic_optimal = cluster_range[np.argmin(aic_scores)]
    sil_optimal = cluster_range[np.argmax(silhouette_scores)]
    # Use silhouette score as primary criterion (best cluster separation)
    optimal_clusters = sil_optimal
   print(f" Cluster estimation results:")
   print(f" BIC optimal: {bic_optimal}")
   print(f" AIC optimal: {aic optimal}")
   print(f" Silhouette optimal: {sil_optimal}")
   print(f" Selected: {optimal clusters}")
   return optimal_clusters, silhouette_scores[sil_optimal - 2]
def fit(self, spike_waveforms, true_num_neurons=None):
    Fit DTW-GMM model to spike waveforms
    Args:
        spike_waveforms: Tensor of shape (n_spikes, window_size, n_channels)
        true_num_neurons: Optional ground truth for comparison
    print(" DTW-GMM Spike Sorting Pipeline")
   print("="*50)
    start_time = time.time()
    # Convert to numpy if needed
    if torch.is_tensor(spike_waveforms):
        spike_waveforms = spike_waveforms.numpy()
   n_spikes, window_size, n_channels = spike_waveforms.shape
```

```
print(f"Processing {n_spikes} spikes ({window_size} samples ×
⇔{n_channels} channels)")
      # Step 1: Select reference spikes for DTW
      print(" Selecting reference spikes...")
      # Use primary channel (max variance) for reference selection
      primary_channel = np.argmax(np.var(spike_waveforms, axis=(0, 1)))
      primary_waveforms = spike_waveforms[:, :, primary_channel]
      self.reference_spikes = self._select_reference_spikes(primary_waveforms)
      print(f" Selected {len(self.reference_spikes)} reference spikes")
      # Step 2: Compute DTW feature vectors
      dtw_features = self._compute_dtw_features_parallel(primary_waveforms)
      # Step 3: Dimensionality reduction and normalization
      print(" Preprocessing DTW features...")
      # Standardize features
      dtw_features_scaled = self.scaler.fit_transform(dtw_features)
      # PCA for further dimensionality reduction
      dtw_features_pca = self.pca.fit_transform(dtw_features_scaled)
      print(f" DTW features: {dtw_features.shape} → PCA: {dtw_features_pca.
⇒shape}")
                 PCA explained variance: {self.pca.explained_variance_ratio_.
      print(f"

sum():.3f}")
       # Step 4: GMM clustering with automatic cluster estimation
      if true_num_neurons is not None:
          # Use ground truth if available
          n_clusters = true_num_neurons
          print(f" Using ground truth: {n_clusters} clusters")
          silhouette_score_final = None
      else:
          # Estimate clusters using overclustering approach
          n clusters, silhouette score final = self.
→_estimate_clusters_gmm_overclustering(
              dtw_features_pca
          )
       # Fit final GMM
      print(f" Fitting final GMM with {n_clusters} components...")
      self.gmm = GaussianMixture(
          n_components=n_clusters,
          covariance_type='full',
```

```
random_state=42,
          max_iter=200
      )
      self.gmm.fit(dtw_features_pca)
      cluster_labels = self.gmm.predict(dtw_features_pca)
      cluster_probs = self.gmm.predict_proba(dtw_features_pca)
      # Calculate final metrics
      if silhouette_score_final is None:
          silhouette_score_final = silhouette_score(dtw_features_pca,__
end_time = time.time()
      # Results summary
      print(f"\n DTW-GMM Clustering Complete!")
      print(f" Processing time: {(end_time - start_time)/60:.2f} minutes")
      print(f" Clusters found: {len(np.unique(cluster labels))}")
      print(f" Silhouette score: {silhouette_score_final:.3f}")
      print(f" Reference spikes: {len(self.reference_spikes)}")
      print(f" DTW features: {dtw features.shape[1]} dimensions")
                PCA features: {dtw_features_pca.shape[1]} dimensions")
      print(f"
      return {
          'cluster_labels': torch.tensor(cluster_labels),
          'cluster_probabilities': torch.tensor(cluster_probs),
          'dtw_features': torch.tensor(dtw_features),
          'pca_features': torch.tensor(dtw_features_pca),
          'n_clusters': n_clusters,
          'silhouette_score': silhouette_score_final,
          'processing_time': end_time - start_time,
          'reference_spikes': torch.tensor(self.reference_spikes)
      }
  def predict(self, new_spike_waveforms):
      Predict cluster labels for new spike waveforms
      if self.gmm is None:
          raise ValueError("Model not fitted. Call fit() first.")
      # Convert to numpy if needed
      if torch.is_tensor(new_spike_waveforms):
          new_spike_waveforms = new_spike_waveforms.numpy()
      # Extract primary channel
```

```
primary_channel = np.argmax(np.var(new_spike_waveforms, axis=(0, 1)))
        primary_waveforms = new_spike_waveforms[:, :, primary_channel]
        # Compute DTW features
        dtw_features = self._compute_dtw_features_parallel(primary_waveforms)
        # Apply same preprocessing
        dtw_features_scaled = self.scaler.transform(dtw_features)
        dtw_features_pca = self.pca.transform(dtw_features_scaled)
        # Predict using trained GMM
        labels = self.gmm.predict(dtw_features_pca)
        probabilities = self.gmm.predict_proba(dtw_features_pca)
       return torch.tensor(labels), torch.tensor(probabilities)
# Integration function for the existing pipeline
def run_dtw_gmm_spike_sorting(spike_waveforms, true_num_neurons=None,
                             sampling_frequency=30000, n_reference_spikes=50):
   Run DTW-GMM spike sorting integrated with the existing pipeline
   Args:
        spike waveforms: Tensor of shape (n spikes, window size, n channels)
        true_num_neurons: Optional ground truth number of neurons
        sampling frequency: Recording sampling frequency
        n_reference_spikes: Number of reference spikes for DTW
    Returns:
       Dictionary with clustering results compatible with existing pipeline
   print("\n DTW-GMM SPIKE SORTING INTEGRATION")
   print("="*60)
   print("Research basis:")
   print("- PMC4749467: DTW for spike classification with temporal tolerance")
   print("- Nature Sci Reports 2019: GMM overclustering approach")
   print("- Hybrid: DTW features + GMM clustering")
   print()
    # Initialize DTW-GMM sorter
   window size = spike waveforms.shape[1]
    sorter = DTWGMMSpikeSorter(
       window size=window size,
       sampling_frequency=sampling_frequency,
       n_reference_spikes=n_reference_spikes,
        dtw_radius=3, # From PMC4749467: small radius for efficiency
        n_jobs=4
```

```
# Fit the model
   results = sorter.fit(spike_waveforms, true_num_neurons=true_num_neurons)
    # Create templates from clusters (compatible with existing pipeline)
   cluster labels = results['cluster labels']
   n_clusters = results['n_clusters']
    # Compute cluster templates
   templates = torch.zeros((n clusters, window size, spike waveforms.shape[2]))
   for cluster_id in range(n_clusters):
       mask = cluster_labels == cluster_id
       if mask.sum() > 0:
            templates[cluster_id] = spike_waveforms[mask].mean(dim=0)
    # Add templates to results for compatibility
   results['templates'] = templates
   results['sorter'] = sorter
   return results
# Example usage integrated with your existing pipeline:
# After extracting spike waveforms in your existing code:
print("\n Running DTW-GMM Spike Sorting...")
# Run DTW-GMM clustering
dtw_gmm_results = run_dtw_gmm_spike_sorting(
   spike_waveforms,
    true_num_neurons=TRUE_NUM_NEURONS,
    sampling_frequency=SAMPLING_FREQUENCY
# Extract results
dtw_cluster_labels = dtw_gmm_results['cluster_labels']
dtw_templates = dtw_gmm_results['templates']
dtw_silhouette = dtw_gmm_results['silhouette_score']
print(f" DTW-GMM Results:")
print(f" Clusters: {dtw qmm results['n clusters']}")
print(f" Silhouette: {dtw_silhouette:.3f}")
# Plot DTW-GMM templates
fiq, axs = plot_waveforms(dtw_templates, name="DTW-GMM template")
fig.suptitle("DTW-GMM Inferred Templates", y=1.05)
plt.show()
```

```
# Compare with ground truth (same evaluation as existing pipeline)
matching_spike_ids = match_true_and_inferred_spikes(
    detected_spike_times, true_spike_times, true_spike_ids
# Create confusion matrix for DTW-GMM
cm_dtw = confusion_matrix(matching_spike_ids, dtw_cluster_labels)
# Plot DTW-GMM confusion matrix
plt.figure(figsize=(10, 8))
plt.imshow(cm_dtw, interpolation="none", cmap=plt.cm.Blues)
plt.colorbar()
plt.xlabel('DTW-GMM Cluster IDs')
plt.ylabel('True Neuron IDs (-1 = False Positive)')
plt.title("DTW-GMM Spike Sorting: Confusion Matrix")
for i in range(cm_dtw.shape[0]):
   for j in range(cm_dtw.shape[1]):
        plt.text(j, i, str(cm_dtw[i, j]), ha='center', va='center',
                color='white' if cm_dtw[i, j] > cm_dtw.max() / 2 else 'black')
plt.tight_layout()
plt.show()
# Calculate DTW-GMM performance
valid_matches = matching_spike_ids >= 0
if valid matches.sum() > 0:
    ari dtw = adjusted rand score(matching spike ids[valid matches],
                                 dtw_cluster_labels[valid_matches])
   print(f" DTW-GMM Performance:")
   print(f" Adjusted Rand Index: {ari_dtw:.3f}")
   print(f" Silhouette Score: {dtw_silhouette:.3f}")
11 11 11
```

[38]: '\n# After extracting spike waveforms in your existing code:\nprint("\n Running DTW-GMM Spike Sorting...")\n\n# Run DTW-GMM clustering\ndtw\_gmm\_results = run\_dtw\_gmm\_spike\_sorting(\n spike\_waveforms, \n true\_num\_neurons=TRUE\_NUM\_NEURONS,\n sampling\_frequency=SAMPLING\_FREQUENCY\n)\n\n# Extract results\ndtw\_cluster\_labels = dtw\_gmm\_results[\'cluster\_labels\']\ndtw\_templates = dtw\_gmm\_results[\'templates\']\ndtw\_silhouette = dtw\_gmm\_results[\'templates\']\n\nprint(f" DTW-GMM Results:")\nprint(f" Clusters: {dtw\_gmm\_results[\'n\_clusters\']}")\nprint(f" Silhouette: {dtw\_silhouette:.3f}")\n\n# Plot DTW-GMM templates\nfig, axs = plot\_waveforms(dtw\_templates, name="DTW-GMM template")\nfig.suptitle("DTW-GMM Inferred Templates", y=1.05)\nplt.show()\n\n# Compare with ground truth (same

```
evaluation as existing pipeline)\nmatching_spike_ids =
match_true_and_inferred_spikes(\n
                                     detected_spike_times, true_spike_times,
true_spike_ids\n)\n\m# Create confusion matrix for DTW-GMM\ncm_dtw =
confusion matrix(matching spike ids, dtw_cluster_labels)\n\n# Plot DTW-GMM
confusion matrix\nplt.figure(figsize=(10, 8))\nplt.imshow(cm_dtw,
interpolation="none", cmap=plt.cm.Blues)\nplt.colorbar()\nplt.xlabel(\'DTW-GMM
Cluster IDs\')\nplt.ylabel(\'True Neuron IDs (-1 = False
Positive)\')\nplt.title("DTW-GMM Spike Sorting: Confusion Matrix")\n\nfor i in
range(cm dtw.shape[0]):\n
                            for j in range(cm_dtw.shape[1]):\n
plt.text(j, i, str(cm_dtw[i, j]), ha=\'center\', va=\'center\',\n
color=\'white\' if cm_dtw[i, j] > cm_dtw.max() / 2 else
\'black\')\n\nplt.tight_layout()\nplt.show()\n\n# Calculate DTW-GMM
performance\nvalid matches = matching spike ids >= 0\nif valid matches.sum() >
        ari_dtw = adjusted_rand score(matching_spike ids[valid matches], \n
0:\n
dtw_cluster_labels[valid_matches])\n
                                       print(f" DTW-GMM Performance:")\n
         Adjusted Rand Index: {ari_dtw:.3f}")\n
print(f"
                                                    print(f"
                                                                Silhouette
Score: {dtw_silhouette:.3f}")\n'
```

```
[41]: print("\n Running DTW-GMM Spike Sorting...")
      # Run DTW-GMM clustering
      dtw gmm results = run dtw gmm spike sorting(
          spike_waveforms_fixed,
          true_num_neurons=TRUE_NUM_NEURONS,
          sampling_frequency=SAMPLING_FREQUENCY
      )
      # Extract results
      dtw_cluster_labels = dtw_gmm_results['cluster_labels']
      dtw_templates = dtw_gmm_results['templates']
      dtw_silhouette = dtw_gmm_results['silhouette_score']
      print(f" DTW-GMM Results:")
      print(f" Clusters: {dtw_gmm_results['n_clusters']}")
      print(f"
                Silhouette: {dtw silhouette:.3f}")
      # Plot DTW-GMM templates
      fig, axs = plot_waveforms(dtw_templates, name="DTW-GMM template")
      fig.suptitle("DTW-GMM Inferred Templates", y=1.05)
      plt.show()
      # Compare with ground truth (same evaluation as existing pipeline)
      matching_spike_ids = match_true_and_inferred_spikes(
          detected_spike_times, true_spike_times, true_spike_ids
      # Create confusion matrix for DTW-GMM
```

```
cm_dtw = confusion_matrix(matching_spike_ids, dtw_cluster_labels)
# Plot DTW-GMM confusion matrix
plt.figure(figsize=(10, 8))
plt.imshow(cm_dtw, interpolation="none", cmap=plt.cm.Blues)
plt.colorbar()
plt.xlabel('DTW-GMM Cluster IDs')
plt.ylabel('True Neuron IDs (-1 = False Positive)')
plt.title("DTW-GMM Spike Sorting: Confusion Matrix")
for i in range(cm dtw.shape[0]):
    for j in range(cm_dtw.shape[1]):
        plt.text(j, i, str(cm_dtw[i, j]), ha='center', va='center',
               color='white' if cm_dtw[i, j] > cm_dtw.max() / 2 else 'black')
plt.tight_layout()
plt.show()
# Calculate DTW-GMM performance
valid_matches = matching_spike_ids >= 0
if valid_matches.sum() > 0:
    ari_dtw = adjusted_rand_score(matching_spike_ids[valid_matches],
                                dtw_cluster_labels[valid_matches])
    print(f" DTW-GMM Performance:")
             Adjusted Rand Index: {ari dtw:.3f}")
    print(f"
    print(f" Silhouette Score: {dtw silhouette:.3f}")
 Running DTW-GMM Spike Sorting...
 DTW-GMM SPIKE SORTING INTEGRATION
______
Research basis:
- PMC4749467: DTW for spike classification with temporal tolerance
- Nature Sci Reports 2019: GMM overclustering approach
- Hybrid: DTW features + GMM clustering
 DTW-GMM Spike Sorting Pipeline
_____
Processing 4708 spikes (416 samples × 32 channels)
 Selecting reference spikes...
 Selected 50 reference spikes
 Computing DTW features: 4708 spikes × 50 references...
DTW features: 100%
                       | 4708/4708 [2:01:36<00:00, 1.55s/it]
 Preprocessing DTW features...
 DTW features: (4708, 50) → PCA: (4708, 9)
  PCA explained variance: 0.952
```

Using ground truth: 8 clusters

Fitting final GMM with 8 components...

DTW-GMM Clustering Complete!

Processing time: 121.97 minutes

Clusters found: 8

Silhouette score: 0.433 Reference spikes: 50

DTW features: 50 dimensions PCA features: 9 dimensions

DTW-GMM Results:

Clusters: 8

Silhouette: 0.433

	DTW	GMM terbal	m	68IM terN€	MASSIM tembl	M	SERIM terRes	NASKIM terbili	M-SIEM 1	DEW	team terbah	<b>te≅</b> M temp
channel 0 [2-score]	100			-r-	ANGEMM tempi							
	100		1			. 			<u></u>	7		L
channel 1 [2-score]	-100										<u></u>	
channel 2 [2-score]	0											
	100	L	]		<u></u>				<u></u>	7	\\ 	\
channel 3 [2-score]	-100 -		1									<u> </u>
channel 4 [2-score]	0								-	-		
channel 5 [2-score]	100											
	100	<u> </u>	]		<u></u>	١.	<u> </u>		<u>L</u>	ᆜ	<u></u>	Ļ.,.,
channel 6 [2-score]	-100		1				·		γ			-r-
channel 7 [2-score]	100	γ					-γ		Lγ			γ
	100	l	]		L	l .	L,		<u></u>	7	L	L
channel 8 [2-score]	-100				<u></u>				_γ	_		
channel 9 [2-score]	0	r	1				·		_γ	_		γ
	100		1									
	-100											
channel 11 [2-score]	0-			-r-					-		<u></u>	
channel 12 [2-score]	100											
	100	L	1			l .			<u></u>	_	닏	Ļ
channel 13 [2-score]	-100										<u></u>	
channel 14 [2-score]	100								_			
	100		]				ļ		<u></u>	<u> </u>		L
	-100											
channel 16 [2-score]	0									-		
	100		1]	—— ——			<b></b>		<u></u>	7		<u> </u>
	-100 -											
channel 18 [2-score]	0	γ							-γ	-		γ
channel 19 [2-score]	100											
	-100 -	L.,'_								_		L.,
channel 20 [2-score]	-100		j						_	_		
channel 21 [2-score]	100											
	100	Ļ	]					$\square$	<u></u>	4	$\square$	<u></u>
channel 22 [2-score]	-100			-v-								
channel 23 [2-score]	100		]									
	100	l, , , , ,				l .	<u> </u>		Ļ.,	7		L
channel 24 [2-score]	-100											
channel 25 [2-score]	0 -		1									
	100					i .	<b></b>		<u></u>	7	$\Box$	<u> </u>
	-100											<u> </u>
channel 27 [2-score]	100		1						v	_		
	100		]									L
	100									]		
channel 29 [2-score]	0 -								-	_		
channel 30 [2-score]	100		]				-29					
	-100 -								<u></u>	_]		<u></u>
channel 31 [2-score]	-100											
	700	-5 0 5 time [ms]		-5 0 5 time [ms]	-5 0 S time [ms]		-5 0 5 time [ms]	-5 0 5 time [ms]	-5 0 time [	ms]	-5 0 5 time [ms]	-5 0 5 time [ms]

```
NameError
                                                 Traceback (most recent call last)
      Cell In[41], line 25
           22 plt.show()
            24 # Compare with ground truth (same evaluation as existing pipeline)
       ---> 25 matching_spike_ids = match_true_and_inferred_spikes(
                   detected_spike_times, true_spike_times, true_spike_ids
           27 )
            29 # Create confusion matrix for DTW-GMM
            30 cm_dtw = confusion_matrix(matching_spike_ids, dtw_cluster_labels)
      NameError: name 'match_true_and_inferred_spikes' is not defined
[47]: # Create SpikeInterface sorting object directly from DTW-GMM results
      detected_spike_times = si_spike_times
      n_spikes = len(dtw_gmm_results['cluster_labels'])
[48]: # Extract DTW-GMM results
      cluster labels = dtw gmm results['cluster labels']
      templates = dtw_gmm_results['templates']
      n_clusters = dtw_gmm_results['n_clusters']
      silhouette_score = dtw_gmm_results['silhouette_score']
      print(f" DTW-GMM Results:")
      print(f" Clusters: {n_clusters}")
                Silhouette: {silhouette_score:.3f}")
      print(f"
      print(f"
                Cluster labels shape: {cluster_labels.shape}")
      # Create dummy spike times if original detection times aren't available
      # This assumes spikes are evenly distributed (for demonstration)
      if 'detected_spike_times' not in locals():
          print(" Creating dummy spike times for demonstration...")
          # Create evenly spaced spike times across the recording
          total_samples = NUM_SAMPLES
          spike_times_dummy = torch.linspace(WINDOW_SIZE, total_samples - WINDOW_SIZE,
                                           len(cluster_labels), dtype=torch.long)
      else:
          # Use actual detected spike times (aligned with waveforms)
          spike_times_dummy = detected_spike_times[:len(cluster_labels)]
      print(f"Spike times shape: {spike_times_dummy.shape}")
```

DTW-GMM Results:

```
Silhouette: 0.433
        Cluster labels shape: torch.Size([4708])
     Spike times shape: torch.Size([4708])
[54]: # Alternative approach: Use ground truth spike times for cleaner evaluation
      def evaluate_dtw_gmm_on_ground_truth_times(dtw_gmm_results, true_spike_times,_u
       →true_spike_ids, traces):
          HHHH
          Evaluate DTW-GMM by applying it to ground truth spike locations
          print(" Evaluating DTW-GMM on ground truth spike times...")
          # Extract waveforms at ground truth locations
          cluster_labels = dtw_gmm_results['cluster_labels']
          sorter = dtw_gmm_results['sorter']
          # Get ground truth waveforms
          gt waveforms = []
          valid_gt_times = []
          valid_gt_ids = []
          for i, t in enumerate(true_spike_times):
              start idx = t - WINDOW SIZE // 2
              end_idx = t + WINDOW_SIZE // 2
              if start_idx >= 0 and end_idx < traces.shape[0]:</pre>
                  waveform = traces[start_idx:end_idx]
                  gt_waveforms.append(waveform)
                  valid_gt_times.append(t)
                  valid_gt_ids.append(true_spike_ids[i])
          gt_waveforms = torch.stack(gt_waveforms)
          valid_gt_times = torch.tensor(valid_gt_times)
          valid_gt_ids = torch.tensor(valid_gt_ids)
          print(f" Extracted {len(gt_waveforms)} ground truth waveforms")
          # Apply DTW-GMM clustering to ground truth waveforms
          gt_predictions, _ = sorter.predict(gt_waveforms)
          # Create SpikeInterface sorting for ground truth evaluation
          sorting_dtw_gt = se.NumpySorting.from_samples_and_labels(
              valid_gt_times.numpy(),
              gt_predictions.numpy(),
              sampling_frequency=SAMPLING_FREQUENCY
          )
```

Clusters: 8

```
# Compare with ground truth
    comparison_gt = sc.compare_sorter_to_ground_truth(
        sorting_true,
       sorting_dtw_gt,
       exhaustive_gt=True,
       match_score=0.1,
       chance_score=0.05
   )
   return comparison_gt, gt_predictions, valid_gt_ids
# Run ground truth evaluation
comparison_gt, gt_pred, gt_true = evaluate_dtw_gmm_on_ground_truth_times(
   dtw_gmm_results, true_spike_times, true_spike_ids, traces
# Visualize ground truth evaluation
fig, ax = plt.subplots(figsize=(10, 8))
sw.plot_confusion_matrix(comparison_gt, ax=ax)
plt.title('DTW-GMM on Ground Truth Spike Times')
plt.show()
```

Evaluating DTW-GMM on ground truth spike times... Extracted 4802 ground truth waveforms Computing DTW features:  $4802 \text{ spikes} \times 50 \text{ references...}$ 

DTW features: 0%| | 19/4802 [02:41<11:19:10, 8.52s/it]

DTW features: 0%| | 11/4802 [05:14<38:00:31, 28.56s/it]

```
KeyboardInterrupt

Traceback (most recent call last)

Cell In[54], line 55

52     return comparison_gt, gt_predictions, valid_gt_ids

54  # Run ground truth evaluation

---> 55 comparison_gt, gt_pred, gt_true = evaluate_dtw_gmm_on_ground_truth_time (

56     dtw_gmm_results, true_spike_times, true_spike_ids, traces

57 )

59  # Visualize ground truth evaluation

60 fig, ax = plt.subplots(figsize=(10, 8))

Cell In[54], line 34, in evaluate_dtw_gmm_on_ground_truth_times(dtw_gmm_results

$\times \true_spike_times, \true_spike_ids, \traces)

31 print(f" Extracted {len(gt_waveforms)} ground \truth waveforms")

33  # Apply DTW-GMM clustering to ground \truth waveforms
```


```
---> 34 gt_predictions, _ = sorter.predict(gt_waveforms)
     36 # Create SpikeInterface sorting for ground truth evaluation
     37 sorting_dtw_gt = se.NumpySorting.from_samples_and_labels(
            valid_gt_times.numpy(),
            gt predictions.numpy(),
     39
     40
            sampling_frequency=SAMPLING_FREQUENCY
     41 )
Cell In[38], line 285, in DTWGMMSpikeSorter.predict(self, new spike waveforms)
    282 primary_waveforms = new_spike_waveforms[:, :, primary_channel]
    284 # Compute DTW features
--> 285 dtw features = self. compute dtw features parallel(primary waveforms)
    287 # Apply same preprocessing
    288 dtw_features_scaled = self.scaler.transform(dtw_features)
Cell In[38], line 106, in DTWGMMSpikeSorter. compute dtw features parallel(self
 ⇔spike_waveforms)
    103
           return dtw_row
    105 # Parallel DTW computation
--> 106 dtw features = Parallel(n jobs=self.n jobs, backend='threading')(
 delayed(compute dtw row)(i) for i in tqdm(range(n spikes), desc="DTW feat res")
    110 return np.array(dtw_features)
File ~/PycharmProjects/STATS320 FinalProject/.venv/lib/python3.12/site-packages
 ⇔joblib/parallel.py:2072, in Parallel._call_(self, iterable)
   2066 # The first item from the output is blank, but it makes the interpreter
   2067 # progress until it enters the Try/Except block of the generator and
   2068 # reaches the first `yield` statement. This starts the asynchronous
   2069 # dispatch of the tasks to the workers.
   2070 next(output)
-> 2072 return output if self.return_generator else list(output)
File ~/PycharmProjects/STATS320 FinalProject/.venv/lib/python3.12/site-packages
 →joblib/parallel.py:1682, in Parallel._get_outputs(self, iterator, pre_dispato_)
   1679
   1681
            with self._backend.retrieval_context():
-> 1682
                yield from self._retrieve()
   1684 except GeneratorExit:
           # The generator has been garbage collected before being fully
   1685
   1686
            # consumed. This aborts the remaining tasks if possible and warn
   1687
            # the user if necessary.
            self._exception = True
   1688
File ~/PycharmProjects/STATS320 FinalProject/.venv/lib/python3.12/site-packages
 →joblib/parallel.py:1800, in Parallel._retrieve(self)
   1789 if self.return_ordered:
```


```
# Case ordered: wait for completion (or error) of the next job
   1790
   1791
            # that have been dispatched and not retrieved yet. If no job
   (...)
   1795
            # control only have to be done on the amount of time the next
            # dispatched job is pending.
   1796
   1797
            if (nb jobs == 0) or (
   1798
                self. jobs[0].get status(timeout=self.timeout) == TASK PENDING
   1799
            ):
-> 1800
                time.sleep(0.01)
                continue
   1801
   1803 elif nb_jobs == 0:
   1804
            # Case unordered: jobs are added to the list of jobs to
   1805
            # retrieve `self._jobs` only once completed or in error, which
   (...)
            # timeouts before any other dispatched job has completed and
   1811
   1812
            # been added to `self._jobs` to be retrieved.
KeyboardInterrupt:
```

```
[53]: | # FAST SOLUTION: Use existing DTW-GMM results without recomputing
      def quick_spikeinterface_evaluation(dtw_gmm_results):
          Quick evaluation using existing DTW results without recomputation
          print(" Quick SpikeInterface evaluation (no DTW recomputation)")
          # Extract existing results
          cluster_labels = dtw_gmm_results['cluster_labels']
          n_clusters = dtw_gmm_results['n_clusters']
          # Create dummy spike times for SpikeInterface (evenly distributed)
          n_spikes = len(cluster_labels)
          spike_times_dummy = torch.linspace(
              WINDOW_SIZE,
              NUM_SAMPLES - WINDOW_SIZE,
              n_spikes,
              dtype=torch.long
          )
          print(f"Using {n_spikes} spikes with {n_clusters} clusters")
          # Create SpikeInterface sorting object
          sorting dtw = se.NumpySorting.from samples and labels(
              spike_times_dummy.numpy(),
              cluster labels.numpy(),
              sampling_frequency=SAMPLING_FREQUENCY
```

```
# Load ground truth
   recording, sorting true = se.read_mearec("data/sim_data/recordings2.h5")
    # Create comparison (this is fast)
    comparison = sc.compare_sorter_to_ground_truth(
       sorting_true,
       sorting_dtw,
       exhaustive_gt=True,
       match_score=0.1,
       chance_score=0.05
   )
   return comparison
# Run the quick evaluation (should take seconds, not minutes)
comparison_quick = quick_spikeinterface_evaluation(dtw_gmm_results)
# Visualize immediately
import spikeinterface.widgets as sw
fig, ax = plt.subplots(figsize=(10, 8))
sw.plot_confusion_matrix(comparison_quick, ax=ax)
plt.title('DTW-GMM Results - Quick Evaluation')
plt.show()
# Get performance metrics
performance = comparison_quick.get_performance()
confusion_matrix = comparison_quick.get_confusion_matrix()
print(" DTW-GMM Quick Evaluation:")
print(f"Confusion Matrix:\n{confusion_matrix}")
if len(performance) > 0 and 'accuracy' in performance.columns:
    accuracy = performance['accuracy'].mean()
   print(f"Average Accuracy: {accuracy:.3f}")
print(f" Evaluation completed in seconds (not minutes)!")
```

Quick SpikeInterface evaluation (no DTW recomputation) Using 4708 spikes with 8 clusters





DTW-GMM Quick Evaluation:

~	_						
Cor	1 🕇 11	21	on	Ma	1.7	יוֹי צ	•

	0.0	1.0	2.0	3.0	4.0	5.0	6.0	7.0	FN
#0	0	0	0	0	0	0	0	0	247
#1	0	0	0	0	0	0	0	0	281
#2	0	0	0	0	0	0	0	0	394
#3	0	0	0	0	0	0	0	0	225
#4	0	0	0	0	0	0	0	0	879
#5	0	0	0	0	0	0	0	0	733
#6	0	0	0	0	0	0	0	0	1122
#7	0	0	0	0	0	0	0	0	924
FP	180	1144	135	134	176	2395	414	130	0
۸		۸	0	000					

Average Accuracy: 0.000

Evaluation completed in seconds (not minutes)!

```
[57]: def robust_spike_detection_fixed(recording, detection_threshold=8.0,
                                     min_isi_ms=3.0, extraction_buffer_ms=0.5):
          Fixed spike detection with proper refractory period and alignment
          Based on search results: "Template matching detected ~85-90% of spikes
          compared to ~70% for the standard fixed threshold method"[7]
          traces = recording.get_traces()
          n_samples, n_channels = traces.shape
          fs = recording.sampling_frequency
          # Convert time parameters to samples
          min_isi_samples = int(min_isi_ms * fs / 1000) # 2ms = 60 samples at 30kHz
          buffer_samples = int(extraction_buffer_ms * fs / 1000) # 0.5ms buffer
          print(f" Fixed spike detection:")
          print(f"
                     Min ISI: {min_isi_ms}ms ({min_isi_samples} samples)")
                     Buffer: {extraction_buffer_ms}ms ({buffer_samples} samples)")
          all_spike_times = []
          for ch in tqdm(range(n_channels), desc="Processing channels"):
              ch_trace = traces[:, ch]
              # Robust noise estimation using MAD
              noise_mad = np.median(np.abs(ch_trace - np.median(ch_trace)))
              noise_std = noise_mad / 0.6745
              threshold = detection_threshold * noise_std
              print(f"
                         Channel {ch}: threshold = {threshold:.2f}")
              # Find all points above threshold
              above_threshold = np.where(np.abs(ch_trace) > threshold)[0]
              if len(above_threshold) == 0:
                  continue
              # FIXED: Proper spike detection with refractory period
              spike times = []
              i = 0
              while i < len(above_threshold):</pre>
                  current_time = above_threshold[i]
                  # Skip if too close to previous spike (refractory period)
```

```
if spike_times and current_time - spike_times[-1] < min_isi_samples:
                i += 1
                continue
            # Find local maximum in a small window around current detection
            search_start = max(0, current_time - buffer_samples)
            search_end = min(n_samples, current_time + buffer_samples + 1)
            search_window = ch_trace[search_start:search_end]
            if len(search window) > 0:
                # Find the actual peak (maximum absolute value)
                abs_window = np.abs(search_window)
                local_max_idx = np.argmax(abs_window)
                actual_spike_time = search_start + local_max_idx
                # Verify this is actually a local maximum
                peak_window_start = max(0, actual_spike_time - 3)
                peak_window_end = min(n_samples, actual_spike_time + 4)
                peak_window = np.abs(ch_trace[peak_window_start:
 →peak_window_end])
                if (len(peak_window) > 0 and
                    np.abs(ch_trace[actual_spike_time]) == np.max(peak_window)):
                    spike_times.append(actual_spike_time)
                    all_spike_times.append(actual_spike_time)
            # Skip ahead past the refractory period
            while i < len(above_threshold) and above_threshold[i] <=__
 current_time + min_isi_samples:
                i += 1
    # Remove duplicates and sort
   all_spike_times = sorted(set(all_spike_times))
   print(f" Fixed detection: {len(all_spike_times)} spikes")
   return torch.tensor(all_spike_times)
# Apply the fixed detection
print(" Applying fixed spike detection...")
fixed_spike_times = robust_spike_detection_fixed(recording_filt,
                                                detection_threshold=4.0, #_
 → Higher threshold
                                                min_isi_ms=1.0) # 2ms_
 ⇔refractory period
```

```
# Re-extract waveforms with fixed detection
fixed_spike waveforms, fixed_valid events = extract_waveforms_robust(
    recording_filt.get_traces(), fixed_spike_times, WINDOW_SIZE
print(f"Fixed extraction: {len(fixed_spike_waveforms)} clean waveforms")
# Plot the first 5 fixed waveforms to verify
if len(fixed spike waveforms) > 0:
    fig, axs = plot_waveforms(fixed_spike_waveforms[:5], name="fixed waveform")
    fig.suptitle("Fixed Spike Detection - First 5 Waveforms", y=1.05)
    plt.show()
 Applying fixed spike detection...
 Fixed spike detection:
  Min ISI: 1.0ms (32 samples)
  Buffer: 0.5ms (16 samples)
                                  | 3/32 [00:00<00:01, 20.61it/s]
Processing channels:
                      9%1
  Channel 0: threshold = 23.00
  Channel 1: threshold = 22.61
  Channel 2: threshold = 22.42
  Channel 3: threshold = 22.40
  Channel 4: threshold = 22.27
Processing channels: 19%|
                                 | 6/32 [00:00<00:01, 23.46it/s]
  Channel 5: threshold = 22.56
  Channel 6: threshold = 23.21
  Channel 7: threshold = 24.20
  Channel 8: threshold = 24.68
Processing channels: 28%|
                            | 9/32 [00:00<00:01, 21.18it/s]
  Channel 9: threshold = 24.00
Processing channels: 38%|
                            | 12/32 [00:00<00:00, 20.72it/s]
  Channel 10: threshold = 23.70
  Channel 11: threshold = 22.75
  Channel 12: threshold = 22.47
  Channel 13: threshold = 22.31
Processing channels: 47% | 15/32 [00:00<00:00, 22.27it/s]
  Channel 14: threshold = 22.19
                          | 18/32 [00:00<00:00, 23.24it/s]
Processing channels: 56%
  Channel 15: threshold = 22.47
  Channel 16: threshold = 22.78
  Channel 17: threshold = 23.06
```

```
Channel 18: threshold = 23.49
   Channel 19: threshold = 23.52
Processing channels: 75% | 24/32 [00:01<00:00, 22.98it/s]
  Channel 20: threshold = 22.95
   Channel 21: threshold = 22.51
  Channel 22: threshold = 22.75
  Channel 23: threshold = 22.31
  Channel 24: threshold = 22.16
Processing channels: 84%
                           | 27/32 [00:01<00:00, 24.12it/s]
  Channel 25: threshold = 22.01
   Channel 26: threshold = 22.11
  Channel 27: threshold = 22.84
  Channel 28: threshold = 22.58
   Channel 29: threshold = 22.58
Processing channels: 94% | 30/32 [00:01<00:00, 24.35it/s]
  Channel 30: threshold = 22.68
Processing channels: 100% | 32/32 [00:01<00:00, 23.11it/s]
   Channel 31: threshold = 22.59
 Fixed detection: 25706 spikes
Extracting waveforms with window size 416...
Skipping spike 0: out of bounds (start=-206, end=210)
Skipping spike 1: out of bounds (start=-201, end=215)
Skipping spike 2: out of bounds (start=-200, end=216)
Skipping spike 3: out of bounds (start=-193, end=223)
Skipping spike 25698: out of bounds (start=1919586, end=1920002)
Skipping spike 25699: out of bounds (start=1919587, end=1920003)
Skipping spike 25700: out of bounds (start=1919588, end=1920004)
Skipping spike 25701: out of bounds (start=1919687, end=1920103)
Skipping spike 25702: out of bounds (start=1919701, end=1920117)
Skipping spike 25703: out of bounds (start=1919702, end=1920118)
Skipping spike 25704: out of bounds (start=1919788, end=1920204)
Skipping spike 25705: out of bounds (start=1919789, end=1920205)
 Extracted 25694 valid waveforms
Fixed extraction: 25694 clean waveforms
```

	fixe	d wavefofi	melli wavefofi	meži wavefofin	œ <b>ž</b> i wavefofin	se≝ waveform 4
channel 0 [2-score]	100	volijkala-javive	physplania	stobjesiger	upolicingum.	att folgospiane
	100					
channel 1 [2-score]	-100 ·		- Internation	- Manada / James	- Innerthead	- Instales factors
channel 2 [z-score]	-100	-commonwealth-alphai	duspysiantyny	- Margal Aurigroupsy	- Martin Andrews	- Harplandy myster
channel 3 [2-score]	100			death of shoot	Seed of bear	
	-100					
channel 4 [2-score]	-100	entitation and the c	-manufallpro	danskipter	-surphilynus	- Amelophilyrous
channel 5 [2-score]	100 ·	unyd yndyfryghaf,	and refra	processor and	manifest from	sering/sens
channel 6 [2-score]	100	Laboration of the laboratory o		audadā a		
	-100					
channel 7 [2-score]	-100	(naneklination			mun frans	
channel 8 [2-score]	100 ·	politica prije		many/ma	- market frequence	
channel 9 [2-score]	100 ·	positiv-orthodox-or	and a second play	- Section of Physics	hereightere	palmont/shops
channel 10 [2-score]	100					
	-100					
channel 11 [2-score]	-100	laid-vestrostilo	new September	- verificatelyheire	- validative	- vitrige-statisfulare-se
channel 12 [2-score]	100 ·	yaphajirmayavasi.	temperous-p-vi	aderioalpi-volun	- sdeeppalp Jungan	sijosaingb/p-vikeni
channel 13 [2-score]	100	unahr.m.				
	-100 100					
channel 14 [2-score]	-100	-manifestration-	i vangagaylada			
channel 15 [2-score]	100	Continuido (sector	evapoly.4/we	ngul/Numan	ngudy.4/wwqqa	updy/freed
channel 16 [2-score]	100		1			
	-100			THEY STATE OF	100,000	Laction 1
channel 17 [2-score]	-100	-resolvent syste	-	-warrant war	-war-wy-are	-
channel 18 [2-score]	100	to-statement	-erromant/su	www.chinage	stomped/subje	Asrageri (fuelle-
	-100 ·		<u></u>		H	
channel 19 [2-score]	-100	aphydayanelun	- Appropriately to	stephen/,heer	- negited, hear	engins/,heer
channel 20 [2-score]	100	pol-l-weighwise	all the second second	- Aprobadji kropina	- April 19 Amples	Aprilospinopino
	-100 100		ļ.,			
channel 21 [2-score]	-100	physion of the same	-edinology-y-lik	Sempleson	Storythisty spilitals	- Marchine Control
channel 22 [2-score]	100	Natural and State of the State	- specially special	-interpretation	- Vicina/Proprint plants	- Guara (registry de
channel 23 c [2-score]	100	and the same of th		de de la	day to the same	in the same
	-100					
channel 24 [2-score]	-100	postajo optija	-	- myspaethings	- majoritoritalismo	- more property deposits
channel 25 [2-score]	100	- Angerest (ma	riorestation	prograph apple	gradiniskusysia.	protection of the second
channel 26 cf [2-score]	100	allegram proposeds.		annya/nam		
	-100 100		Ļ	-		
channel 27 [2-score]	-100	antichistor		- mary prome	-von/h/hom	- Applement
channel 28 [2-score]	100 ·	aletantist factorian	-stratefylpu	patrici farma	Naple Lanne	Application of the same of the
channel 29 [2-score]	100	enderly it regions	deservence	- 6-4-5-00//freein	- www.	
	-100 ·		Ļ	<del></del>	<u></u>	
channel 30 [2-score]	0 -100	lalpolpopography	desiration of the	45~		
channel 31 [2-score]	100	Anddellermily	innsheeds/lyt	nshamilitelijisini	ushamidi(dythic	- malanan deli di filipi di ma
d S	-100	-5 0 5 time [ms]	-5 0 5 time [ms]	-5 0 5 time [ms]	-5 0 5 time [ms]	-5 0 5 time [ms]

```
[59]: # Use MEArec's known spike times for validation
      recording, sorting true = se.read mearec("data/sim_data/recordings2.h5")
      # Get ground truth spike times
      true_spike_times = []
      for unit_id in sorting_true.unit_ids:
          unit_spikes = sorting_true.get_unit_spike_train(unit_id)
          true_spike_times.extend(unit_spikes)
      true spike times = sorted(true spike times)
      print(f"Ground truth: {len(true_spike_times)} spikes")
      # Extract waveforms at KNOWN spike locations
      def extract_at_known_locations(traces, spike_times, window_size):
          waveforms = []
          for spike_time in spike_times[:100]: # First 100 for testing
              start_idx = int(spike_time) - window_size // 2
              end_idx = int(spike_time) + window_size // 2
              # Fix: Check against traces.shape[0] for number of samples
              if start idx >= 0 and end idx < traces.shape[0]:</pre>
                  # Fix: Extract all channels using [:, :] indexing
                  waveform = traces[start idx:end idx, :]
                  waveforms.append(waveform)
          # Fix: Convert to torch tensor for plot_waveforms compatibility
          if len(waveforms) > 0:
              return torch.tensor(np.array(waveforms), dtype=torch.float32)
          else:
              return torch.tensor([], dtype=torch.float32)
      # Extract waveforms at ground truth locations
      gt_waveforms = extract_at_known_locations(
          recording_filt.get_traces(),
          true_spike_times,
          WINDOW SIZE
      )
      print(f" Extracted {len(gt_waveforms)} ground truth waveforms")
      print(f"Ground truth waveform shape: {gt_waveforms.shape}")
      # Plot these to see what real spikes should look like
      if len(gt_waveforms) > 0:
          fig, axs = plot_waveforms(gt_waveforms[:5], name="ground truth")
```

Ground truth: 4805 spikes
Extracted 99 ground truth waveforms
Ground truth waveform shape: torch.Size([99, 416, 32])
DTW features: 4%| | 199/4802 [57:51<22:18:14, 17.44s/it]

		ound truth	ground truth	ground truths	round truth	ground truth
channel 0 [2-score]	0 -	productive production of the second	- species plan	hospitalend	. Applier of the process	- Arranament
	100 -			<u></u>	<b>!</b>	<b>!</b>
channel 1 [2-score]	-100 -	anderston/an	- Miles diles (Virtuals	- product/product	, i s a constantino	- Andreadour
channel 2 [2-score]	100 -	el-cooperate transfer to	Harabi Andy Nedley	helpf-harfy/redistro	PMENDERPY	*************
	100 -			-	$\vdash$	<u></u>
channel 3 [2-score]	-100 -	umaphasa kayatin Ang	dianity (type)	- insplayor provide	- the section of	- Market plants
channel 4 [2-score]	100 -	-andrew (militare	-	whitem	-2014/04/2014/4/4/10	-
	-100 -			-	<u></u>	<u></u>
channel 5 [2-score]	100 - 0 - -100 -	mended/an	- sundystyrus		- Palarantanggar	-
	100 -		mar. 14.			
channel 6 [z-score]	0 - -100 -	mentaly	teriso (fun)	anner (Arrie	- managharian hyper	- Landerson,
channel 7 [2-score]	0 -				namayr	merhan
	100 -			+		<b>1</b> ——
channel 8 [2-score]	0 - -100 -	managha	-market / harme	and from	minimal pool (*	Lung
channel 9 [2-score]	100 -	and the state of t	Undergriji Stepher		why	MM
	100 -				$\square$	<u> </u>
channel 10 [2-score]	0 -	ytherprepations			- mariprotroma	
	100 -					
channel 11 [2-score]	-100 -	nume op regulates	- medical deliner	e-sp-uselyter-cor	Antonic (Specialism	- age-served Ag
channel 12 [2-score]	100 -	tourposessive of	- ninesembrane	- processing-viscous	-yer-busing-at	-
	100 -			+	$\vdash$	<u></u>
channel 13 [2-score]	-100 -	terendigi-sphijzhel	and give the frequency	- magazaphi, indones	energi alijih jeryi	waysherme
channel 14 [2-score]	100 -	t-co-grang-yluria	n-grymy-ylvinoyd	-grang-ylaskophys	- vight plight plant part	Spellperspelsous,
	100 -				$\square$	
channel 15 [2-score]	0 -	unqua <sub>r</sub> ajuus	- AND A TOWARD	ant of most	- whombindum, to	-
	100 -					
	0 - -100 -	n-spikulalivie	man Agrana	- And Andrewson		and and a
channel 17 [2-score]	100 -	menye	- maryen	-HT-WANE	- manager of the	marapay
	100 -			1	<del></del>	1,
channel 18 [2-score]	-100 -	-enventer/su	stamped frame	-	who are the fire	- hondulands
channel 19 [2-score]	100 -	-April property of the	mingstoned, freen	Mphwlytonia		mym
	-100 -				H	<u> </u>
channel 20 [2-score]	0 -	ng.having.kraj	- April of London	ender@Longlands	maryman	- mylm
	100 -					
	-100 -		Threshill and philipse		- marchato,	honh
channel 22 [2-score]	0 -	market and the second second	- share (my Hydra-in	- and the physical and	-	
	100 -		<del>   </del>   <del>  </del>	+	<del></del>	1
channel 23 [2-score]	-100 -	mindered type gates	-dank-farjanjon	JAA-'ge'ylayed	frankrakanja.	-
		riceningueside	- marking a septiment of the septiment o	printerstations	-decel-deceler	Lympunus
					$\vdash$	<u></u>
channel 25 [2-score]	0 -	rigrasinations	- prospresión región		- edwards-drives	Inter-designary
	100 -					
	-100 -	-deriendu (estr	Eventus/Herver		posterjenjevi	
channel 27 [2-score]	100 -		-may popular	mple//man	- jngendjerane-ytt	present-fatheds
				-	$\vdash$	<del></del>
channel 28 [2-score]	-100 -	www.	neightenn	saylallamod	sprije nikovalute	-dunal, it justs
		AND CONTRACTOR	anning finan	~04/48/5/\$nation	Strantoviva	abanny mpany
€2.	100 -				H	<u></u>
Score	0 -	allegerijs-halleler-M <sub>e</sub> n <sup>t</sup> hyse		48	- revisionally of	manyagan
	100 -					
	0 - -100 -	-5 0 S	-S 0 S time [ms]	-S 0 S	-5 0 5 time [ms]	-S 0 S
		-5 0 5 time [ms]	time [ms]	time [ms]	time [ms]	-5 0 5 time [ms]

Ground truth statistics:

Amplitude range: 34.0 - 142.4

Mean amplitude: 101.3 Std amplitude: 34.0

[]:[